

# Pray in memory from linked content

Open a Rosary from the memorial with the linked person already set as the intention.

6 min · Available when Devotions are enabled for the community · 2026-09-21

## Before you begin

✓ Open a published memorial entry.


✓ Check that Devotions are available.

1

STEP 1 OF 3 · 2 MIN

## Open the memorial

Open the community memorial and choose a published person. The prayer link appears only when Devotions are available and the memorial entry is public.




**Fig. 1** The starting point is an entry already visible in the memorial.

2

STEP 2 OF 3 · 2 MIN

## Select Pray in their memory

Open the entry and select Pray in their memory. The command resolves the published person and opens the Rosary; it does not light a candle or begin prayer yet.



**Fig. 2** The link is available because the entry is published in the memorial.


3

STEP 3 OF 3 · 2 MIN

## Review the intention

On the preparation screen, check that the intention reads In memory of Giovanni Verdi. You can change the intention mode before beginning.

The figure stops before Begin: no prayer has been started or completed.



**Fig. 3** The linked intention is ready, but prayer has not begun.

## What you have achieved

- ✓ You can verify the linked name before beginning the Rosary.

## Troubleshooting

### The link does not appear



The entry must be published in the memorial and Devotions must be enabled.